

Hj-6A Thermostatic Magnetic Hotplate Stirrer/Electric Mixer, Hj-6, Hj-^B



Basic Info

- Model NO.: HJ-6A
- Trademark: JY/OEM
- Origin: Changzhou, China
- Type: Magnetic Stirrer
- Specification: 1 pc in 1 box
- HS Code: 8479820090

Product Description

Multi-Unit Magnetic Hotplate Stirrer/Electric Mixer

1. 4 units, thermostatic, LED display speed and temp, easy to use, 4 units are available
2. ISO13485/14001/9001, SGS Audited
3. Since 1984

Magnetic Hotplate Stirrer, an ideal essential device for laboratories, is widely applied in colleges and universities, environment protecting, researching, sanitation, antilipidemic, petroleum, metallurgy, chemical and medical units.

Specifications of HJ-6A Multi-unit Thermostatic Magnetic Stirrer Hotplate	
Power Supply	AC 220V/110V 50Hz/60Hz
Speed	0-2400rpm, LED display
Motor power	25W*6
Heating power	150W*6
Timer	120m
Temperature control range: room temp to 100° C, LED display	

The instrument is equipped with temperature setting and measuring the speed step-less adjustable in order to have and the temperature, speed, digital display, is intelligent temperature magnetic stirrer apparatus. The instrument is widely used among major colleges, environmental protection, scientific research, health and epidemic prevention, chemical, medical and other units in the laboratories for stirring and mixing the liquid solution.