

JJ-3A Digital Display Four-in-One Laboratory Electric Stirrer



Basic Info

- Model NO.: JJ-3A
- Structure: Desktop
- Certification: ISO, SGS
- Type: Electric Stirrer
- Specification: 1 pc in 1 box

- Customized: Customized
- Material: Steel
- Application: School, Hospital, Lab
- Trademark: JY
- Origin: Changzhou, China

Product Description

Electric Stirrer

1. Hi-power, multitask, four unit, timer, separate speed, digital display
2. ISO13485/14001/9001 Certified, SGS Audited
3. Since 1984

The DC digital-display constant speed power-driven stirrer has sound design and novel structure. Adopting DC brushless motor and step-less speed governing with a stabilizing system, it has advantages of digital speed display, high power, low noise and easy operation, and is a basic equipment for mixing chemical reactions with various kinds of liquid phases in laboratories in all universities, colleges, scientific and research institutes and chemical plants.

Specifications of JJ-3A Digital Display Four-in-one Laboratory Electric Stirrer	
Power Supply	AC 220/110V 50/60Hz
Power output	4*100W
Speed	0-1800rpm, step-less regulation
Timer	0-120min
Speed control	Separate control, digital display speed