JJ-4 Four-in-One Synchronous Electric Stirrer/Laboratory Stirrer



Basic Info

Model NO.: JJ-4Stirrer: Six-Unit

Specification: 1 pc in 1 boxHS Code: 8479820090

• Type: Electric Stirrer

• Trademark: JY

• Origin: Changzhou, China

Product Description

Electric Stirrer

1. Hi-power, multitask, four unit, timer, synchronous speed

2. ISO13485/14001/9001 Certified, SGS Audited

3. Since 1984

The DC digital-display constant speed power-driven stirrer has sound design and novel structure. Adopting DC brushless motor and step-less speed governing with the stabilizing system, it has advantages of digital speed display, high power, low noise and easy operation, and is a basic equipment for mixing chemical reactions with various kinds of liquid phases in laboratories in all universities, colleges, scientific and research institutes and chemical plants.

Power Supply	AC 220/110V 50/60Hz
Motor power	200W
Speed	0-3000rpm
Timer	0-120min
Digital display	for speed