

JJ-5 Temperature Control Electric Stirrer/Mixer



Basic Info

- Model NO.: JJ-5
- Structure: Desktop
- Certification: ISO, SGS
- Type: Electric Stirrer
- Specification: 1 pc in 1 box

- Customized: Customized
- Material: Steel
- Application: School, Hospital, Lab
- Trademark: JY
- Origin: Changzhou, China

Product Description

Laboratory Stirrer

1. LED display temp control, timer
2. ISO13485/14001/9001, SGS audited factory
3. Since 1984 The DC digital-display constant speed power-driven stirrer has sound design and novel structure. Adopting DC brushless motor and step-less speed governing with a stabilizing system, it has advantages of digital speed display, high power, low noise and easy operation, and is an basic equipment for mixing chemical reactions with various kinds of liquid phases in laboratories in all universities, colleges, scientific and research institutes and chemical plants.

Specifications of JJ-5 Temperature Control Electric Stirrer/Mixer	
Power Supply	AC 220/110V 50/60Hz
Motor power	100W
Speed	0-3000rpm
Timer	0-120min
Digital display	for temperature
Temp control	RT-100celsius
Synchronous drive, step-less speed regulation, easy to operate.	