

JJ-6 DC Motor Constant-Speed Laboratory Stirrer/Mixer



Basic Info

- Model NO.: JJ-6
- Structure: Desktop
- Certification: ISO, SGS
- Type: Electric Stirrer
- Transport Package: Size: 46L*34W*14h (Cm) G.W/N.W: 10kg/9kg
- Origin: Changzhou, China

- Customized: Customized
- Material: Steel
- Application: School, Hospital, Lab
- Trademark: JY
- Specification: 1 pc in 1 box
- HS Code: 847982009

Product Description

Laboratory Stirrer

1. LED, Timer, DC motor, Constant-speed
2. ISO13485/14001/9001, SGS audited factory
3. Since 1984The DC digital-display constant speed power-driven stirrer has sound design and novel structure. Adopting DC brushless motor and step-less speed governing with a stabilizing system, it has advantages of digital speed display, high power, low noise and easy operation, and is an basic equipment for mixing chemical reactions with various kinds of liquid phases in laboratories in all universities, colleges, scientific and research institutes and chemical plants.

Specifications of JJ-6 DC Motor Constant-speed Laboratory Stirrer	
Power Supply	AC 220/110V 50/60Hz
Motor power	100W/200W/300W
Speed	0-2600rpm
Timer	0-120 with hold function

Stepless speed regulation, easy to operate, high quality.