


JY-JX5L Vertical gel electrophoresis cell



Technical specifications

Glass plate size (W×L) : 216×220 mm

Gel Size (W×L) : 186×205 mm

Gel thickness: 1.0 (mm)

Gel amount: 1~2 (pieces)

Sample volume: (1.0mm) 25、40、52wells

Buffer Volume: ~3500 (ml)

Dimension (L×W×H) : 280×140×250 (mm)

Net weight: 5.6 (kg)

Characteristic feature

- The installation of the glass plate in just 15 seconds, and does not need the knob, fast and convenient.
- With the gel making device modular, you can make multi slab gel at the same time.
- Safety cover open button design, convenient opening the cover.
- Abundant buffer not only guarantees cooling effects but also keeps the pH Value stable during the whole experiment process.
- Auto-switch-off when the lid is removed.