

JY-SCZ2+ Mini-Protean Vertical Electrophoresis Cell



Technical specifications

Glass plate size (W×L) : 100×100 mm

Gel Size (W×L) : 82×88 mm

Gel thickness: 0.75 、 1.0、 1.5 (mm)

Gel amount: 1~2 (pieces)

Sample volume: (0.75mm) 11、 15 wells ; (1.0mm) 11、 15wells ; (1.5mm) 11、 15wells

Buffer Volume: ~750 (ml)

Dimension (L×W×H) : 150×120×115 (mm)

Net weight: 1.5 (kg)

Characteristic feature

- The installation of the glass plate in just 15 seconds, and does not need the knob, fast and convenient.
- With a gel casting in the original position, no need to move glasses from gel making to electrophoresis running. It is more convenient to observe gel preparing through both sides of the glasses.
- High transparent polycarbonate material injection molding, impact resistance, high-temperature resistance, corrosion resistance.
- Safety cover open button design, convenient opening the cover.
- With the background color, easy to be observed in the sample adding an electrophoresis process.
- Abundant buffer not only guarantees cooling effects but also keeps the pH Value stable during the whole experiment process.
- Auto-switch-off when the lid is removed.