

JY-SP7 Cellulose Acetate Membrane Electrophoresis



Technical specifications

Bridge Size (W×L): 70×250mm×2; 90×250mm×2.

Buffer Volume: ~1000ml.

Dimension (L×W×H) : 350×260×65 (mm)

Net weight: 1.6 (kg)

Characteristic feature

- The JY-SP7 is updated from the JY-SP7A with a double-bridge style design;
- With a cooling device, which can make the electrophoresis process under constant temperature;
- You can adjust the position of the adjusting rod according to the medium length, to meet the different electrophoretic distance.
- Unfolded electrodes make the repair and parts replacement easy and quick.
- Auto-switch-off when the lid is removed.