

# JY-SPCT Horizontal Electrophoresis



## Technical specifications

Gel Size ( W×L ) : 60×60mm; 120×60mm; 60×120mm 120×120mm.

Sample volume: 3, 6,8,11,13, 18 ,25.

Buffer Volume: ~550ml.

Dimension ( L×W×H: 300×160×75 ( mm )

Net weight: 1.1 ( kg )

## Characteristic feature

- Abundant buffer not only guarantees cooling effects but also keeps the pH Value stable during the whole experiment process.
- Unfolded electrodes make the repair and parts replacement easy and quick.
- With the sample background color design on the tray, convenient to adding the sample.
- The transparent top lid effectively prevents liquid in the cell from volatilizing and avoids electricity leakage.
- Auto-switch-off when the lid is removed.