


# JY-SPFT Horizontal Electrophoresis Cell



## Technical specifications

Gel Size ( W×L ) : 130×130mm; 78×100mm; 48×75mm.

Sample volume: 3, 5, 10, 15, 13, 26.

Buffer Volume: ~600ml.

Dimension ( L×W×H : 305×152×80 ( mm )

Net weight: 1.7 ( kg )

## Characteristic feature

- Abundant buffer not only guarantees cooling effects but also keeps the pH Value stable during the whole experiment process.
- Unfolded electrodes make the repair and parts replacement easy and quick.
- With the fluorescent ruler on the tray, can be observed directly into the ultraviolet.
- With a special gel casting, preventing the gel leakage, the operation is simple and convenient.
- The transparent top lid effectively prevents liquid in the cell from volatilizing and avoids electricity leakage.
- Auto-switch-off when the lid is removed.