

JY-SX3 2D gel Electrophoresis system



Technical specifications

ID Tube Gel Electrophoresis:

length of Glass tube: 120mm

Gel diameter: 1.0, 2.0(mm)

2D Vertical Electrophoresis

Gel size (W×L) : 143×139mm

Gel thickness: 1.0 (mm), 2.0 (mm)

Gel Amount: 1~2 (pieces)

Sample volume: (1.0mm) 21 、 29wells ; (2.0mm) 21 、 29wells

Buffer Volume: ~ 2500 (ml)

Volume (L×W×H) : 460×230×400 (mm)

Weight: 13.2kg