

JY-ZY7 Western Blotting System



Technical specifications

Trans area (W×L) : 216×110 mm

Buffer Volume: ~ 1800 (ml)

Dimension (L×W×H): 290×170×155 (mm)

Net weight: 1.7 (kg)

Characteristic feature

- With the color marking on the splint, to ensure put the gel right into the inner core.
- Unfolded electrodes make the repair and parts replacement easy and quick.
- The transparent top lid effectively prevents liquid in the cell from volatilizing and avoids electricity leakage.
- Auto-switch-off when the lid is removed.
- Not only can be used as a device independent, and can be used as a module with the buffer tank of JY-JX5 and JY-CZ-B.