

JY300HC Universal Power Supply

Technical specifications

Type of Output: Constant-Voltage,
Current or Power

Increment: 1V, 1mA, 1W

Output Range: 5–300V, 1–1500mA,
1–450W.

Timer Range: 1min~99hr, 59min

Display: LCD display with backlight

Output jacks: Four sets of output
jacks



Characteristic feature

- Molding shell, touch keys, dual core microprocessor intelligent control;
- Indicate the preset value and the actual output value at the same time;
- It can store 10 electrophoresis methods;
- With automatic memory function;
- With standard, timing operation function;
- With constant voltage, constant current, constant power, misoperation, fault intelligent prompt functions;
- Automatic detection of no-load, over-load, short circuit, rapid resistance change, ground leak and system overheating.