

# KB-3-D Spiral Mixer



This is a small mixer that can be used in the cold storage house. It is suitable for mixing of blood collection. The mixing movement of this mixer is three-dimensional and does not produce bubbles. It is the ideal choice for a sample with a small dose. The concave pad of the KB – 3 – D mixer can firmly fix tubes. Applicable tubes are 60 x 1.5ml, 30 x 3 ml, 21 x 7 ml and 17 x 15 ml.

Rotate speed:	6 rpm
Motion / obliqueness:	3 – D / fixed 20 degrees
Maximum load:	0.8KG
Working temperature:	4 to 65 degrees Celsius
Power supply:	230 V / 50 HZ
Overall dimensions:	241 x 172 x 172 mm
Weight:	1.2 kg