

## KJMR-III Multi-Mixer



### Description:

: KJMR-III Multi-Mixer has the capability to perform three modes of mixing, all of which are programmable. It is very simple to operate within seconds. It can mix variable substance such as blood sample, urine sample, injectant, blood plasma etc.

### Mixing method:

1. Rotation: Rotate 360 , clockwise, anticlockwise or clockwise and anticlockwise alternately.
2. Vibration: Shake up and down.
3. Vortex: Swirl the samples according to the mode of vortex.

### Specification:

Voltage : 220v $\pm$ 10%50/60Hz

Speed: 5-30r/m

Working plate1: 1.5ml-15ml tube in 22 pcs

Working plate2: 5 0ml tube in 6 pcs, 1.5ml-15ml tube in 4 pcs as spare parts

Moving angle: 15°-90°(per 15° as a shift) Vortex angle: 1°-5°

Maximum working time: 12 hours

External Dimension: 380\*182\*190mm