

L-150A Series Cold Light Illumination



- More stable power supply with CE standard electrical equipment and PC control circuit.
- With delay startup design to protect eyes and electrical equipment caused by sudden brightness.
- Enhancing the security with skillful power off device.
- Cooling the illumination better with advanced wind loop design.
- Stable structure with heavy deadweight of illuminator.
- Longer working life and lower noisy.

Specification :

		L-150A-1	L-150A-2	L-150A-O
Power Supplier	Switch Power 110V/ 220V			
	T-H Lamp with Reflector 20V/ 150W			
	Brightness Adjustable			
	Multi-coated IR Cutoff Filter for Heat Insulation	●	●	●
	Fan Cooling			
	Size 260×120×150mm			
	Weight 3.3 Kg			
Single Light Guide	Rigid Flectional Tube, Length 500mm	●		
Dual Light Guide	Rigid Flectional Tube, Length 500mm		●	
Ring Light Guide	Flexible Tube, Length 1000mm, Adapter Ring Size Φ50mm/ Φ60mm			●
Kohler Illumination Attachment	Attached on Single or Dual Light Guide to adjust and obtain Even Illumination	○	○	

Note: ● Standard Outfit, ○ Optional