

LX-600 Mini-Centrifuge



Special new design centrifuge rotor with 0.2 ml x 8 x 4 PCR.

Standard configuration micro – centrifuge rotor with 1.5 / 2 ml x 8.

Standard configuration rotor adapter with 0.5 ml x 8 and 0.2 ml x 8.

Special design click mechanism that is easy to insert and switch the centrifuge rotor, no need for any other tools.

Fashionable and small by design and a space savor.

Design philosophy is strong and long – lasting with many years of experience in manufacturing the mini-centrifuges, guarantee durability and high quality

Design to fit many different kinds of micro – centrifuge tubes, upon closing, the rotor starts and the rotor will stop once the cover is opened, easy to operate.

Configuration: 1.5 / 2 ml x 8 micro – centrifuge tube rotor

0.2 ml x 8 x 4 PCR rotor

0.5 x 8 adapter and 0.2 x 8 adapter

Power:	220V 50HZ
Frequency:	8000 rpm
Centrifuge force:	3000g
Dimensions:	155×175×125mm
Speed control:	Auto – constant Speed
Drive mold:	Direct drive