

Triple Beam Balance

The MB2610 and MB2610 Bare Single-pan balance of unequal arm type with 3beams. The maximum capacity is 610 grams without extra attachment weights or 2610 with the supplied attachment weights. The readability is 0.1g and the knife is made of high quality alloy steel and the bearing plane is of a gate, which can improve accuracy of use and lengthen life of operation.

The Mb2610B has a counter weight for tare of up to 225g. This product is widely used for weighing in industrial enterprises (such as chemical, medical and machine-building enterprises), agricultural applications and schools.

Specification:

Features:

- 1) Three beam direct reading.
- 2) Magnetic damping system Outer packing.



Model	MB2610	MB2610B
Capacity	2610g	
Tare range		0~225gonMB2610Bonly
Readability	0.1g	
Beam Calibration	FrontBeam10*0.1g	
	2"Beam500*100g	
	3" Beam100*10g	
Damping	Magnetic	
Pan Size	150mm(6")diameter Stainlesssteel	
Overall size(w*h*d)	470*160*150mm (18.5**6.2**6")	
Packing dimensions(L*W*H)	515*385*180(4pcs)	
Gross weight	14KGS(4pcs)	