



MC-1

MC-1 Multifunctional CCD Camera

MC-1/MD-1 Multifunctional Camera is equipped with color photography device and adjustable photography objective, can capture picture, text,near or far scene. Matching relevant eyepiece adapter, it can be applied in biology microscope or stereo microscope, can display clear and balanced dynamic image on TV or monitor, it's the ideal equipment for multimedia teaching, meeting, network communication and production demo etc.

Spec:

Model	MC-1A	MC-1B
Color sensor	1/3" CCD Horizontal resolution: 380 TV lines	1/3" CCD Horizontal resolution: 520 TV lines
Video output	1Vp-p. 75Ω. NTSC/PAL(Optional)	
White balance	Automatic	
Input power	220V/110V(Optional)	
Video cable	Standard Video cable, length:2.9m	

Model	MD-1A	MD-1B
Color sensor	1/2" CMOS	
Resolution	1.3 Mega pixels (1280X1024) Dynamic in real time USB2.0 Port.	2.0 Mega pixels (1600X1200) Dynamic in real time USB2.0 Port.
White balance	Automatic/Manual	Automatic
System Requirements	Windows 2000/XP/Vista	Windows 2000/XP
USB Cable	USB cable, length:1.8M	