

## MC-1080



- **Standard Configuration**

Model	MC-1080
Eyepiece	Wide field WF10X ( $\mu$ 18mm)
Objective	Achromatic objectives 10X/0.25
Video monocular	Built-in clear color 1/3" CMOS sensor, horizontal resolution: 380TV lines, White Balance: Automatic, Color.
Video cable	Standard Video cable, length:1.8m
Stage	Fixed stage with slide moving holder moving range:70mmX30mm
Focus mechanism	Coaxial coarse/fine focus minimum division of fine focusing: 4 $\mu$ m
Light source	LED base without rechargeable batteries, input: 220V/110V (optional)

- **Optional accessories**

Name	Sort/Technique parameter	NO.
Eyepiece	Wide field WF16X ( $\mu$ 11mm)	1051016
Objective	Achromatic objectives 4X/0.10	2010204
	Achromatic objectives 20X/0.40	2010120

	Achromatic objectives 40X/0.65	2010140
Nosepiece	Triple (Frontward ball bearing outer locating)	030001

**Diagram**

