


# MJX-150 Laboratory Mould Incubator, Mould Cultivation Cabinet

## Basic Info

Model NO.: MJX-150

- Structure: Floor Type
- Certification: ISO, SGS
- Type: Incubator
- Specification: 1pc in 1 box
  
- Customized: Customized
- Material: Stainless Steel
- Application: School, Hospital, Lab
- Trademark: JY/OEM
- Origin: Changzhou, China



## Product Description

150L Mold Incubator

1. 150L, LED, 5000LX, 5-50celsius, UV sterilization, reliable
2. ISO13485/14001/9001 Certified, SGS Audited Factory
3. Since 1984

Temperature, time control, touch type keyboard set the adjustment; Equipped with illumination device in chamber, it is convenient for observation through wide angle of view heat preservation vacuum toughened glass; (Mold incubator is equipped with UV light) floodlight is standard disposition for biochemical/shaker incubator (it can be exchanged to UV light according to the requirement of user for sterilization); User can decide the number of shelves in chamber and adjust them optionally.

<b>Specifications of MJX-150 Mold Incubator</b>	
Power Supply	AC 220V/110V 50Hz/60Hz
Volume	150L
Temperature Range	5-50celsius
Temperature fluctuation	± 0.5celsius
Humidity Range	45-90%
Humidity fluctuation	± 5% RH
Illuminance	5000LX
Digital Display	Yes
Ultraviolet disinfection	Yes
Humidity mode	Ultrasonic, even
Ventilation mode	Air duct
Widely applicable for lab, biology, chemistry, physical, hospital, food and environmental protection fields.	