

MJX-250 Laboratory Thermostatic Mould Cultivation Cabinet

Basic Info

- Model NO.: MJX-250
- Structure: Floor Type
- Certification: ISO, SGS
- Type: Incubator
- Specification: 1pc in 1 box

- Customized: Customized
- Material: Stainless Steel
- Application: School, Hospital, Lab
- Trademark: Jinyi
- Origin: Changzhou, China



Product Description

250L Mold Incubator

1. 250L, LED, 5000LX, 5-50celsius, UV sterilization, reliable
2. ISO13485/14001/9001 Certified, SGS Audited Factory
3. Since 1984

Specifications of MJX-250 Mold Incubator	
Power Supply	AC 220V/110V 50Hz/60Hz
Volume	250L
Temperature Range	5-50celsius
Temperature fluctuation	± 0.5celsius
Humidity Range	45-90%
Humidity fluctuation	± 5% RH
Illuminance	5000LX
Digital Display	Yes
Ultraviolet disinfection	Yes
Humidity mode	Ultrasonic, even
Ventilation mode	Air duct
Widely applicable for lab, biology, chemistry, physical, hospital, food and environmental protection fields.	