



Microscope Cover Slips-Economy

Description	Thickness			Packaging
	Nr.1	Nr.1.5	Nr.2	
8*8mm	✓	✓	✓	10x100pcs/ vacuum pack
10*10mm	✓	✓	✓	
12*12mm	✓	✓	✓	
18*18mm	✓	✓	✓	
20*20mm	✓	✓	✓	
22*22mm	✓	✓	✓	
24*24mm	✓	✓	✓	
24*32mm	✓	✓	✓	
24*40mm	✓	✓	✓	
24*50mm	✓	✓	✓	
24*60mm	✓	✓	✓	
Every 10 boxes wrapped in one Aluminum foil with vacuum				

Microscope Cover Slips-Standard


Description	Thickness			Packing
	Nr.1	Nr.1.5	Nr.2	
18*18mm	✓	✓	✓	200pcs/box
20*20mm	✓	✓	✓	200pcs/box
22*22mm	✓	✓	✓	200pcs/box
24*24mm	✓	✓	✓	200pcs/box
24*32mm	✓	✓	✓	100pcs/box
24*40mm	✓	✓	✓	100pcs/box
24*50mm	✓	✓	✓	100pcs/box
24*60mm	✓	✓	✓	100pcs/box
Individually Wrapped In Vacuum Aluminum Foil.				

Microscope Cover Slips-Circular

Description	Thickness			Packing
	Nr.1	Nr.1.5	Nr.2	
1 9mm	✓	✓	✓	100pcs/ box
1 10mm	✓	✓	✓	
1 12mm	✓	✓	✓	
1 13mm	✓	✓	✓	
1 15mm	✓	✓	✓	
1 16mm	✓	✓	✓	
1 18mm	✓	✓	✓	
1 19mm	✓	✓	✓	
1 20mm	✓	✓	✓	
1 22mm	✓	✓	✓	
1 24mm	✓	✓	✓	
1 25mm	✓	✓	✓	

Cover Glass

Nr.1	Nr.1.5	Nr.2
0.13-0.16mm	0.16-0.19mm	0.19-0.23mm



**Hemocytometer ART.No.1280**

**Counting Chamber Cat No.1103**

Including one set of Counting Chamber, one set of blood diluting pipettes for red corpuscles, one set of blood diluting pipettes for white corpuscles, one set of blood cover galss.

Counting Chamber Cat No.1103

blood diluting pipettes for red corpuscles

blood diluting pipettes for white corpuscles



blood diluting pipettes for red corpuscles

blood diluting pipettes for white corpuscles