

Mt-1 Temperature Adjustable Laboratory Sand Bath



Basic Info

- Model NO.: MT-1
- Structure: Desktop
- Certification: ISO, SGS
- Type: Sand Bath
- Trademark: OEM/Unbrand
- Origin: China Mainland

- Customized: Customized
- Material: Steel
- Application: School, Hospital, Lab
- Metal Body: Thermostatic
- Specification: CE, ISO13485/9001/14001

Product Description

1. Novel design cold-rolled sheet, durable
2. CE, ISO13485/9001/14001
3. Since 1984A sand bath is most commonly used in conjunction with a hot plate or heating mantle. A beaker is filled with sand or metal pellets (called shot) and is placed on the plate or mantle. The reaction vessel is then partially covered by sand or pellets. The sand or shot then conducts the heat from the plate to all sides of the reaction vessel.

Specifications of MT-1 Temperature-adjustable Laboratory Sand Bath	
Power Supply:	AC 220V/110V 50Hz/60Hz
Digital Display:	Yes
Heating power:	2000W
Working capacity:	30L*40W*10D (cm)