

# N-510 Series Multi-Viewing Microscope



- Multi-viewing Microscope allows people observation at the same time.
- With the features of Infinite Optical System, Effective Illumination, LED Pointer and images Coherence, it is widely used in Clinic, Scientific Research and Teaching Demonstration.

## Specification:

		N-510	N-306	N-204	N-204F
Optical System	Infinite Optical System	●	●	●	●
Viewing Head	Compensation Free Trinocular Head, Inclined at 30°, 360° Rotatable, Interpapillary Distance 48-75mm	1PC	1PC	1PC	1PC
	Compensation Free Binocular Head, Inclined at 30°, 360° Rotatable, Interpapillary Distance 48-75mm	4PCS	2PCS	1PC	1PC
Eyepiece	Extra Wide Field Eyepiece EW10×/ 20	10PCS	6PCS	4PCS	4PCS
Nosepiece	Inverse Quintuple Nosepiece	●	●	●	●
Objective	Infinite Plan Achromatic Objective 4×、10×、40×、100×	●	●	●	●
Condenser	Swing Condenser N.A.0.9/ 0.25	●	●	●	●
Focusing	Coaxial Coarse & Fine Focus Adjustment, Fine Division 0.001mm	●	●	●	●
Stage	Double Layers Mechanical Stage 185×142mm, Moving Range 75×55mm	●	●	●	●
Kohler Illumination	External Illumination, Halogen Lamp 24V/ 100W	●	●	●	●
Pointer	Green LED Pointer , Brightness Adjustable	●	●	●	●
	Dual Color Pointer , Brightness Adjustable	○	○	○	○
Photo Attachment		○	○	○	○
Video Adapter	C Mount 1×	○	○	○	○
	C Mount 0.5×	○	○	○	○

Note:● Standard Outfit, ○ Optional