

PH-012 Portable pH Meter



PH-012 Portable pH meter

Large-screen, high resolution and accuracy, reading convenient, fast and accurate. Measuring accurate, simple, rugged and durable, Economical and practical, fully in line with a wide range of industries need to test. The manual 2 point through trimmer to ensure the accuracy of measurement data. BNC with the composite electrode and 9 V battery folded, around the world to facilitate the replacement of batteries and electrodes. In short, a measurement a wide range of high precision, small size, the use of convenient, cheap and good products.

Measuring range	0.00~14.00pH
Accuracy	±0.05PH
Resolution	0.01PH
Operating Temperature	0°C~50°C
Automatic Temperature Compensation	0°C~50°C
Power Supply	1x9V
Dimensions	128mm×68mm×25mm
Calibration	Manual 2 point, through trimmer