

PH-013 Portable pH/mV/Temperature Meter



KL-013 Portable pH/mV/Temperature Meter

KL-013 is a portable pH meter designed to be rugged and easy-to-use, yet be reliable and practical. It is ideal for education and for measurements in the field. Simple pH calibration is performed manually through 2 trimmers. Measurements are automatically compensated for variations in temperature. It shows the reading together with the current mode (pH, mV or °C).

Model	PH-013	
Measuring range	PH	0.00~14.00PH
	mV	-1999~1999mV
	Temp	0~100°C
Accuracy	PH	±0.01pH
	mV	±0.1%F·S±1 digit
	Temp	±0.4°C
Resolution	PH	0.01pH
	mV	1mV

	Temp	0.1°C
Environment	Temperature:0°C~50°C Humidity:≤95%	
Display	4-digital LCD	
Power Supply	9V	
Input Impedance	10 ¹² Ω	
Dimensions	180mm×83mm×46mm	
Weight	290g	