

PH-013M Portable pH Meter



PH-013M Portable pH/mV/Temperature Meter

Use the composite electrode, but also for improving the electrode and after-sales service, so that the life of this equipment has greatly extended, Application of advanced electronic technology, high-quality components, with more than 10 years of production experience, and fully ensure the equipment is in power, measurements very steady, measured the scope and application of wide and high accuracy, easy to operate. KELILONG to the new technology, the lowest price and high quality public service users.

Measuring range	PH	0.00~14.00PH
	mV	0~±1999mV
	°C	0~100°C
Resolution	PH	0.01PH
	mV	1mV
	°C	0.1°C
Accuracy	PH	±0.01pH
	mV	±0.1% F·S±1digit
	°C	±0.4°C

Display	4-digital LCD
Power Supply	1x9V
Environment	Temperature:0°C~50°C Humidity:≤95%
Input Impedance	1012 Ohm
Dimensions	123mm×72mm×33mm
Weight	238g