

PH-8813 Portable pH Meter



PH-8813 Portable pH and Temperature meter

Spec:

Measuring range	PH	0.00~14.00PH
	°C	0~100°C
Resolution	PH	0.01PH
	°C	0.1°C
Accuracy	PH	±0.02PH
	°C	±0.5°C
Display	PH	4-digital LCD
	°C	3-digital LCD
Automatic Temperature Compensation	0°C~50°C	
Power Supply	1x9V	
Operating Temperature	0°C~50°C	
Calibration	2 points with auto buffer recognition	
Dimensions	100mm×146mm×48mm	