

PH-98 Lab. pH/ORP/Temperature Meter



PH-98 Lab. pH/ORP/Temperature Meter

Descriptions:

This instrument and built-in microprocessor chips, good looking, feature rich, easy to use, with the following notable features:

1. High accuracy 0.001
2. new designs, large screen displays
3. User-friendly interface, stability indicating identity
4. Three standard, two-point or single-point automatic calibration
5. Automatic temperature compensation, or manual temperature compensation of
6. Wide range of PH and redox potential measurements

Spec:

Range	pH	-2.000 ~ 16.000 pH
	mV	-1999.9 ~ 1999.9mV
	Temp	0 ~ 100 ° C
Resolution	pH	0.001 pH

	mV	0.1 mV
	Temp	0.1 ° C
Accuracy	pH	± 0.001pH
	mV	± 0.1 mV (0 ~ 550 mV) ; ± 0.2 mV (550 ~ 1999mV)
	Temp	± 0.1 ° C
Automatic Temperature Compensation	0 ~ 50 ° C	
pH Calibration	1 or 2 points with auto buffer recognition	
Input impedance	10 12 Ohm	
Power Supply	12 Vdc	
Environment	Temperature:0 °C ~50 °C Humidity: ≤ 95 %	
Dimensions	150 x 110 x 50 mm	