

# QB-128 Rotational Mixer



It can mix different pipe racks or two of the same pipe racks at one time.  
It has a large capacity with simultaneous mixture and the rotation mixture of small samples.  
A simple fixture devices can be quickly installed or removed.  
Rotational type speed adjustment is simple and practical. It can arbitrarily set the mixing requirement  
of 10 – 80 rev / min.

Power supply:	220V
Capacity:	40W
Applicable temperature range:	4 – 60 degrees Celsius
Mixed angles:	Freely adjustable
Rotational velocity:	10 – 80 turns / min
Rotary speed control:	Step – less speed regulation
Maximum mixing quantity:	1200 ml
Minimum mixing quantity:	1.5 ml
Overall dimensions:	290×220×360mm
Instrument configuration:	one principle machine    two pipe support ( can select two)