

# QB-228 Rolling Incubator



Certified by  
ISO9001 International Quality Control Standard.  
Mixing angle is adjustable with horizontal position for small – amplitude mixing and vertical position for full mixing.  
Mainframe has an anti – spray platform with 3 optional accessories of tube racks applicable to sample tubes in different specs.  
Suspension made to handle organisms or blood samples needing air exchange.  
Digital rotating speed display with timer of 0 – 120 min / continuous with precise control.  
Optimized experimental conditions with high repeat ability.

Power supply:	220V 50HZ
Rotating speed:	10 – 80 rpm
Rotating speed control:	Digital
Timing:	0 – 120 min / continuous
Exterior dimensions:	290×220×360mm