

QB-528 Mini Incubator



Adjustable rotational Angle, horizontal position for slightly mix, vertical position to fully mix.
Speed adjusting knob can be set to 6-80 rpm/min as required.
Two tube holders are suitable for different sample tubes used at the same time.
It can be used for gas exchange regarding biological samples and blood samples.

Power supply: 110-230V

Capacity: 30 W

Temperature range: 4°C-60°C

Mixing angle: 0-90 degrees

Rotational velocity: 6 – 80 rpm / min

Rotate speed control: step-less speed regulation

Overall dimensions: 230×150×230mm

Standard configuration : 50ml*2/1.5ml*4+15ml*2/1.5ml*6 tube holder
