

## QB-600 High Speed Oscillation Mixer



This is our company's latest mixer that was improved with reference to Japanese technology. It will quickly mix and oscillate with multiple containers. The functions of manual or automatic mixing can be achieved by exchanging different fittings. The test tubes in the form of 50 ml, 15 ml, and 1.5 ml will be mixed efficiently bring a top-class mixing oscillator for your experiments.

|                     |                          |
|---------------------|--------------------------|
| Oscillation mode:   | Gyration                 |
| Oscillation speed:  | 300 – 2800 rpm           |
| Gyration amplitude: | 3mm                      |
| Capacity:           | 35 W                     |
| Timing:             | 0 – 120 min / continuous |
| Weight:             | 6.5KG                    |