

QB-8002 Micro-Plate Mini Shaker



New design and elegant, this mini shaker is made with 96 specific microplates.

It is also made for long – lasting.

The motor will run continually without overheating.

ABS plastic molding that is classic and made to last. It is also environmentally friendly and made without any spray paint.

Various voltage available and works between 80 – 240 Voltage.

The platform is designed to fit one or two microplates.

Power:	80/240V 50/60HZ
Capacity:	20 W
Rotation radius:	15 mm
Frequency:	0 – 1500 rpm
Speed:	Step-less speed regulation
Oscillation radius:	2 mm
Platform type:	QB-8002/2 micro-plate
Overall dimensions:	260 x 150 x 110 mm
Weight:	3 kg