

QB-9001 Microporous Quick Shaker



Breaking through convention with first-class originality and updated structure. Streamline mainframes formed by enhanced type engineering plastics, up to 4 to 6 mini w-plate quick shaking, non-skid natural rubber shaking plate, and timing and continuity all together bring you an optimum platform for experiments.

Power:	220V
Capacity:	30W
Speed:	Step-less speed regulation
Frequency:	0 – 1350 rpm
Timing:	0 – 120 min / continuous
Loading:	Four 96-pore mini plate
Overall dimensions:	350×250×150mm