

QL-902 Vortex Mixer



It has a ultra small design with new appearance. The QL – 902 is molded by ABS engineering plastics without any spray of paint. The body is acid, alkalis, and corrosion resistant. It provides point oscillation and continuous function with the natural wear – resistant rubber bowl, which will accelerate the experiment. The cast iron stable base will not move or rotate. It is suitable for quick mixing of different sample cells with 0.5 ml to 50 ml.

Blending speed:	2800 turns / min
Working method:	Vibrating and continuous
Power: V	220
Capacity:	30 W
Overall dimensions:	150 x 105 x 115 mm
Weight:	3 kg