

QL-861

Power: 220V50Hz Capy: 40W Operation Mode: Touching, Continuous Rotary Speed: 3000RPM Plate: Ø55mm Rubber

9-06



Rotation speed: 3000rpm Operation mode: continuous, step less speed regulation Plate:Ø55mmrubberCapy:40wPo wer:220v50hz



QL-901VortexStirrer

Model	QL-901
Speed	2800r/min
Power	220V40W
WorkingMode	Continuous,Touch
WorkingPlatform	¢ 60mmWearablerubberhead
Measure	160*120*170