

# RY-2 MELTING POINT TESTER



## RY-2 MELTING POINT TESTER

Melting point is the temperature of something turning to liquid from solid. Testing it is the main method to detect some characters such as purity etc.

RY series are suitable for testing Melting points of drug, spice, and dye etc.

### Applicable standard

Corporate standard(Thickness tester) Q/12XQ0185-2010

### Features:

Test temperature with a standard mercury thermometer.

Adjust heating voltage by keystroke.

Display voltage with a digital voltmeter.

Observe Melting sample by a magnifier.

Test three samples at once.

### Specifications:

Range of melting point: (50~300)°C

Accuracy:  $\pm 1.0^{\circ}\text{C}$

Usage temperature: (15~35)°C

Rate: (1.0、2.0、4.0、6.0)°C/min

Power of the heater: 100W

Power: 220V/ 50Hz/ 100W

Dimension: (190\*120\*200)mm<sup>3</sup>