

## SME-116M Biological Microscope

Specification:

Specification		Model
		SME-116M
Viewing Head	Sliding Binocular Head, Inclined at 45 °	●
Eyepiece	WF10x	●
Objective	Achromatic 4x, 10x, 40x(s), 100x(s, oil)	●
Stage	Double layers mechanical stage, Stage size:104*95mm, Moving Range:50mm*15mm	●
Condenser	N.A.1.25 Abbe Condenser with iris diaphragm and filter	●
Focusing	Coaxial coarse & fine focusing adjustment, Fine focusing scale value 0.002mm	●
Light Source	Reflected rechargeable LED light, brightness adjustable	●
	Transmitted rechargeable LED light, brightness adjustable	●
Packing	Foam + Carton packed	●
Optional Accessories	Triple nosepiece WF16X, WF20X eyepiece; 20X, 60X (s) objective; Monocular Head; Trinocular Head Incandescent Lamp 220V/20W or 110V/20W Wooden Packed or Aluminium Case packed	



SME-116M