

SME-F Series



Specification				
		SME-F	SME-F1	SME-F1-100
Viewing Head	Monocular Inclined at 45°, 360° Rotatable	●	●	●
Eyepiece	WF10×/18mm	●	●	●
Objective	Achromatic 4×, 10×, 40×(s)	●	●	
	Achromatic 4×, 10×, 40×(s), 100×(s) oil			●
Stage	Plain stage with slide clips Stage size: 110mm×110mm	●		
	Double layers mechanical stage Stage size: 110mm×126mm Moving range: 70mm×30mm		●	●
Condenser	Single lens with Disc diaphragm	●		
	N.A.1.25 Abbe Condenser with iris diaphragm and filter		●	●
Focusing	Coaxial coarse & fine focusing adjustment Fine focusing scale value 0.002mm	●	●	●
Light Source	Incandescent lamp 220V/20W or 110V/20W	●	●	●
Optinal Accessory	Eyepiece: WF16×, WF20×, Objective: 20×, 60×(s) Halogen bulb 6V/20W, AC 220V/110V LED illumination			