

## SME-F4 Series



Specification		Model		
		SME-F4	SME-F4D	SME-F4DT
Viewing Head	Monocular Inclined at 45°, 360° Rotatable	●		
	Compensation Free Binocular Head, Inclined at 30° (55mm-75mm)		●	
	Compensation Free Trinocular Head, Inclined at 30° (55mm-75mm)			●
Eyepiece	WF10×/18mm	●	●	●
Objective	Achromatic 4×, 10×, 40×(s), 100×(s) oil	●	●	●
Stage	Stage size: 120mm×125mm Mechanical stage Moving range: 60mm×30mm	●	●	●
Condenser	N.A.1.25 Abbe Condenser with iris diaphragm and filter	●	●	●
Focusing	Separated coarse & fine focusing adjustment	●	●	●
Light Source	Halogen bulb 6V/20W, AC 220V/110V Adjustable light source	●	●	●
Optional Accessory	Eyepiece: WF16×, WF20×, P16× Objective: 20×, 60×(s) Semi-plan or Plan-achromatic objectives Reflective Mirror, Demonstration (Teaching) Head CCD Adapter Compensation Free Binocular & Trinocular Head (48mm-75mm) LED illumination Fluorescent lamp 220V/5W or 110V/5W			