

SME-F6 Series



Specification		Model		
		SME-F6	SME-F6B	SME-F6B-100
Viewing Head	Monocular Inclined at 45°, 360° Rotatable	●	●	●
Eyepiece	WF10X/18mm	●	●	●
Objective	Achromatic 4X, 10X, 40X(s)	●	●	
	Achromatic 4X, 10X, 40X(s), 100X(s) oil			●
Stage	Plain stage with slide clips Stage size: 110mmX110mm	●		
	Double layers mechanical stage Stage size: 110mmX126mm Moving range: 70mmX30mm		●	●
Condenser	N.A.0.65 Condenser with disc diaphragm	●		
	N.A.1.25 Abbe Condenser with iris diaphragm and filter		●	●
Focusing	Coaxial coarse & fine focusing adjustment Fine focusing scale value 0.002mm	●	●	●
Light Source	Incandescent lamp 220V/20W or 110V/20W	●		
	Halogen bulb 6V/20W, AC220V/110V		●	●
Optimal Accessory	Eyepiece: WF16X, WF20X, Objective: 20X, 60X(s) LED illumination			