

SME-F6BN



Specification		Model		
		SME-F6A	SME-F6BN (SME-F6D)	SME-F6E
Viewing Head	Monocular Inclined at 45°, 360°Rotatable	Y		
	Compensation Free Binocular Head,Inclined at 30°(55mm-75mm)		Y	
	Compensation Free Trinocular Head,Inclined at 30°(55mm-75mm)			Y
Eyepiece	WF10x/18mm		Y	
Objective	Achromatic 4x, 10x, 40x(s), 100x(s, oil)		Y	
Stage	Double layers mechanical stage, Stage size:125*130 mm, Moving Range:70*30 mm		Y	
Condenser	N.A.1.25 Abbe Condenser with iris diaphragm and filter		Y	
Focusing	Coaxial coarse & fine focusing adjustment, Fine focusing scale value 0.002mm, Focusing Range 30mm		Y	
Light Source	Halogen bulb 6V/20W, AC 220V/110V			
Optical Accessory	Eyepiece:WF16x, WF20x, P16x			
	Objective :20x, 60x(s)			
	Reflective Mirror, Demonstration (Teching) Head			
	CCD Adapter			
	Compensation Free Binocular & Trinocular Head(48mm-75mm)			
	LED illumination			