

SME-F9A



SME-F9A



SME-F9B



SME-F9T

Specification		Model		
		SME-F9A	SME-F9B	SME-F9T
Viewing Head	Monocular head, inclined at 45 °, 360° rotatable	●		
	Compensation free binocular head, inclined at 30°(48mm- 75mm)		●	
	Compensation free binocular head, inclined at 30°(48mm- 75mm)			●
Eyepiece	WF10×/18mm	●	●	●
Objective	Achromatic 4×、 10×、 40×(s)、 100×(s) Oil	●	●	●
Stage	Double layer mechanical stage	●	●	●
	Stage size: 126mm×110mm	●	●	●
	Moving Range: 70mm×30mm	●	●	●
Condenser	N.A.1.25 Abbe condenser with iris diaphragm & filter holder	●	●	●
Focusing	Coaxial coarse & fine focusing adjustment, fine focusing scale value 0.002mm	●	●	●
Light Source	LED1W, 220V/110V	●	●	●
Optional Accessory	Eyepiece: WF16×、 WF20×、 P16× Objective: 20× (s)、 60×(s) LED light Semi-plan objective: 4×、 10 ×、 20×、 40×(s)、 60×(s)、 100×(s) oil Plan objective: 4×、 10×、 20×、 40×(s)、 60×(s)、 100×(s) oil Dual viewing head Darkfield and phase contrast attachment 1.3Mega、 3.0Mega、 5.0Mega、 10.0Mega COMS digital camera eyepiece Photography attachment and CCD C-mount: 0.5X、 0.57X、 0.75X			