

SME-F9K



SME-F9K



SME-F9M



SME-F9MB

Specification		Model			
		SME-F9	SME-F9K	SME-F9M	SME-F9MB
Viewing Head	Monocular head,inclined at 45 °, 360° rotatable	•	•	•	•
Eyepiece	WF10 X 18mm	•	•	•	•
Objective	Achromatic: 4X , 10X , 40X(S)	•	•		
	Achromatic : 4X , 10X , 40X(S), 100X (S) Oil			•	•
Stage	Plain stage with slide clips,Stage size: 110mm X 120mm	•	•		
	Mechanical stage Stage size: 110mm x 120mm Moving range: 60mm x 30mm			•	
	Double layers mechanical stage Stage size: 110mm x 126mm Moving range: 70mm x 30mm				•
Condenser	0.65 Condenser with disc diaphragm	•			
	0.65 Condenser with disc diaphragm		•		
	N. A.1.25 Abbe Condenser with iris diaphragm and filter holder			•	•
Focusing	Coaxial coarse & fine focusing adjustment Fine focusing scale value 0.002mm	•	•	•	•
Light Source	LED1W,220V/110V	•	•	•	•
Optional Accessory	Eyepiece: WF16x,WF20x,P16x Objective: 20X、 60X(S) Semi-plan Objective: 4X,10X,20X(S),40X(S),60X(S),100X(S)Oil Plan Objective: 4X,10X,20X(S),40X(S),60X(S),100X(S)Oil Dual Viewing Head 1.3Mega、 3.0Mega、 5.0Mega、 10.0Mega pixels COMS Digital camera eyepiece Photography attachment and CCD C-mount 0.5X、 0.57X、 0.75X				