

TD4N Table top low speed centrifuge



Feature:

- 1. The Frame is 3 tiers protection steel jacket, and with the stainless-steel chamber, inside protective steel lagging which is simple-compact, light-weight and low noise.
- 2. Adopt advanced CPU control system realizing microprocessor control, it can control rotate speed, and relative RCF, digital(LED) display which indicates the current speed, time and RCF.
- 3. Brushless DC motor, no powder pollution, unnecessary to maintain, electric lid lock.
- 4. With many kinds of rotor for your choice.
- 5. Over speed and imbalance protection, the centrifuge will stop automatically when the above situation happened, with short spin function.
- 6. Emergency locks under the bottom of the centrifuge which help you open the lid when power off.
- 7. It is consistent with CE security certification.

Specification:

Max speed	4000rpm	Max RCF	2950×g
Max volume	4×100ml	Power supply	AC 220V 50HZ 2A
Timer	0~99min	Noise	≤55dBA
Dimension (D x W x H)	528×370×280mm	Net weight	28 KG
Speed accuracy	±20rpm	Package	carton box or wooden box

Optional Rotors:

No	Rotors No	Max speed r/min	Volume(ml)	Max RCF(g)
No31277	Swing out rotor	4000	4×4×15ml	2950
No31277	Swing out rotor	4000	4×6×10ml	2580
No31277	Swing out rotor	4000	4×100ml	2820
No31277	Swing out rotor	4000	4×50ml	2820
No30611	Angle rotor	4000	6x50ml	1820