

# TM-1 TRANSPARENCY TESTER



## TM-1 TRANSPARENCY TESTER

Transparency tester is requisite for detecting transparency of gelation.

### Applicable standard

National standard (Food additive Gelatin) GB6783-94

National standard (Hard gelatin capsule for medicine) GB13731-92

Corporate standard (Thickness tester) Q/12XQ0193-2010

### Features:

Thoroughly eliminating the influence to results by environment and experimenter's visual error.

Standard test tube and same capacity of the sample.

High stable light source system, equal brightness.

High precise photo electric transducers for testing exactly.

Auto-testing, transparency being conversion accounted automatically.

Automation: auto-test, auto-diagnose, auto-alarm.

### Specifications:

Test range: (20~500)mm

Test accuracy:  $\pm 5$ mm

Volume: 400ml

Power: 220V/ 50Hz/ 80W

Dimension: (360×240×120) mm<sup>3</sup>