


TM-2 TRANSPARENCY TESTER



TM-2 TRANSPARENCY TESTER

Transparency tester is requisite for detecting transparency of gelation.

Applicable standard

National standard (Food additive Gelatin) GB6783-94

National standard (Hard gelatin capsule for medicine) GB13731-92

Corporate standard (Thickness tester) Q/12XQ0193-2010

Features:

Thoroughly eliminating the influence to results by the environment.

The sample being conveniently taken out or put in by using a standard test tube.

Airtight light source system.

Specifications:

Test range: (20~500)mm

Test accuracy: ± 5 mm

Power: 220V/ 50Hz/ 80W

Dimension: (600×200×630) mm³