



Optional accessories

Vortex-5 Vortex Shaker

Power:220V Capy:50W

Rotation Speed: 2800RPM

Working: Continuous, Point touching,

Speed Regulation

Plate: Exchangeable Between Bowl Type

and Flat Type Overall Dimension:170*120*170mm

Feature: Powerful vortex shaker of point touching, continuity and speed regulation

Mixer mode of touching and vortex: It maybe controlled through speed regulation to mix from low-speed to high-speed with proper modes to

choose as required.

Multifunction: Both bowl type shaking plate and flat type shaking is

provided form annual and automatic mixing modes of different test tubes and

containers;

Stable operation: Enough weight and integrated metal shell allow various

shaking to have a stable operation platform.

VORTEX-6

The Vortex-6 shaker upgraded from VORTEX-5 is uniquely equipped with light-operated device with no need of frequent power supply switching, totally different from original operation mode and having really materialized the forward-minded service state-of-the-art.

Power: 220V Capy:50W

Rotation Speed:2800RPM

Working: Continuous, Point touching, Speed Regulation Plate: Exchangeable Between Bowl Type and Flat Type Overall

Dimension: 170*120*170mm



